

Flügelhorn
1 in B

Bayrische Bläserweisen

Josef Lang jun.

Hinterberger Weis'

Musical notation for the first piece, 'Hinterberger Weis''. It consists of two staves. The first staff starts with a treble clef, a key signature of two flats (B-flat and E-flat), and a 3/4 time signature. The music begins with a *mf* dynamic and includes a *f* dynamic later. The second staff starts at measure 11 and continues the melody with a *f* dynamic.

Musical notation for the second piece, 'Im Kiefernwald'. It consists of three staves. The first staff starts with a treble clef, a key signature of one sharp (F#), and a 3/4 time signature. The music begins with a *f* dynamic. The second staff starts at measure 9 and continues the melody with a *f* dynamic. The third staff starts at measure 19 and continues the melody with a *mf* dynamic.

Copyright (c) 2008 by Edition Rinner / ER-061

Blau-Weis' - Schön is' im Bayernland

Musical score for 'Blau-Weis' - Schön is' im Bayernland'. The score is written in treble clef with a key signature of one flat (B-flat) and a time signature of 3/4. It consists of three staves. The first staff begins with a measure rest of 2 measures, followed by a series of eighth notes. The second staff starts at measure 9 with the instruction 'flott' and features accents (^) over several notes. The third staff starts at measure 17 and includes a first ending bracket and a second ending marked '2.'. Dynamics include 'f' (forte) and 'mf' (mezzo-forte).

Ausnahms-Weis' - Geh ma ins Gebirg

Musical score for 'Ausnahms-Weis' - Geh ma ins Gebirg'. The score is written in treble clef with a key signature of three flats (B-flat, E-flat, A-flat) and a time signature of 2/4. It consists of three staves. The first staff begins with a measure rest of 3 measures, followed by eighth notes. The second staff starts at measure 9 and includes dynamics 'mf' and 'f'. The third staff starts at measure 17 and includes time signature changes to 2/4 and 3/4, along with dynamics 'mf' and 'f'.

Copyright (c) 2008 by Edition Rinner / ER-061